



Certificate of Analysis

Catalog Number: TSL-314-1
Lot Number: 171130-087
Toxin: T-2 Toxin
Quantity: 1.0 mL
Concentration or
Concentration after reconstitution: 100.0µg/mL
Reconstitution if applicable: N/A
Solvent: Acetonitrile
Storage: 8°C or lower
Date of Manufacturing: November 30, 2017
Date of Expiration: 1year from certificate date (see below).

Purity: ≥ 98%
Traceability: All components used in production of this standard can be tracked back to their source for each batch produced, allowing full product traceability.
Testing Method: LC MS/MS
Standard Uncertainty (µg/mL): ±2.24
Product Specification (µg/mL): 90.0-110.0
Results (µg/mL): 101.2

This is to certify that T-2 Toxin standard Lot number 171130-087 passes the required specifications and is released for sale.

Craig Humphrey – Chemistry Production Manager
March 28, 2018