



# Certificate of Analysis

Catalog Number: TSL-314  
Lot Number: 160824-286  
Toxin: T-2 Toxin  
Quantity: 5.0 ml  
Concentration or  
Concentration after reconstitution: 100.0µg/ml of T-2 Toxin  
Reconstitution if applicable: N/A  
Solvent: Acetonitrile  
Storage: 8°C or lower  
Date of Manufacturing August 24, 2016  
Date of Expiration: 1 year from certificate date (see below).

---

Purity: ≥ 98%  
Traceability: All components used in production of this standard can be tracked back to their source for each batch produced, allowing full product traceability.  
Testing Method: LC MS/MS  
Standard Uncertainty (µg/ml): T-2 Toxin: ±2.24  
Product Specification: T-2 Toxin: (90.0µg/ml-110.0µg/ml)  
Results (µg/ml): T-2 Toxin: 102.7µg/ml

This is to certify that T-2 Toxin standard Lot number 160824-286 passes the required specifications and is released for sale.

Carrie Maune – Quality Manager  
October 12, 2016

Trilogy Analytical Laboratory, Inc.  
870 Vossbrink Drive Washington, MO 63090  
Phone: 636-239-1521 Fax: 636-239-1531 [www.trilogylab.com](http://www.trilogylab.com)

REVISION HISTORY

<b>REVISED DATE</b>	<b>REVISED BY</b>	<b>CHANGES MADE</b>
5/11/10	DANIELLE GANT	ORIGINAL DOCUMENT
6/24/10	Danielle Gant	Added page number and code
1/13/11	Danielle Grimes	Corrected Code Number
9/26/2011	Danielle Grimes	Added paginated numbers

Trilogy Analytical Laboratory, Inc.  
870 Vossbrink Drive Washington, MO 63090  
Phone: 636-239-1521 Fax: 636-239-1531 [www.trilogylab.com](http://www.trilogylab.com)