



Certificate of Analysis

Catalog Number: TSL-317
Lot Number: 180612- 220
Toxin: Deoxynivalenol
Quantity: 10.0 mL
Concentration or Concentration after reconstitution: 100.0µg/mL
Reconstitution if applicable: N/A
Solvent: Methanol
Storage: 8°C or lower
Date of Manufacturing: June 12, 2018
Date of Expiration: 1year from certificate date (see below).

Purity: ≥ 98%
Traceability: All components used in production of this standard can be tracked back to their source for each batch produced, allowing full product traceability.
Testing Method: UV Spectroscopy
Testing Equipment Conditions: HP 8452A Diode Array Spectrophotometer
Solvent: Methanol
Wavelength: 219 nm
E Value: 6208
Standard Uncertainty (µg/mL): ±2.24
Product Specification (µg/mL): 90-110
Results (µg/mL): 98.3

This is to certify that DON standard Lot number 180612- 220 passes the required specifications and is released for sale.

Craig Humphrey – Chemistry Production Manager
August 8, 2018